

**FOR IMMEDIATE RELEASE**

May 12, 2006



Contact: Dana Peterson, Public Relations, (503) 922-5597

**Portland, Ore** ~ Banfield, The Pet Hospital® introduces its patented Anesthesia Safety Pressure Relief Valve. The valve's unique pop-off design helps ensure pets will no longer face anesthetic complications or death due to a safety pressure relief valve being left shut. The valve is designed to stay open, and can only shut when it is forcefully pressed down. It is universally adaptable to virtually any anesthesia machine.

Includes:

- Safety Pressure Relief Valve
- (1) Stop Sign Instruction Card with accessories
- (1) Installation Instructions
- (1) Machine and Vaporizer Leak Test Card  
New instructions for doing a high pressure leak test. Please ensure that your surgical team is aware of how the new equipment works and how to perform a high pressure leak test BEFORE beginning a surgery
- (1) Pop-off Valve Functional Settings Card

The Safety Pressure Relief Valve is available exclusively through MWI and is also available on LEI Medical anesthesia machines. To place an order, contact MWI at 800-824-3703, reference MWI part #028362, [www.mwivet.com](http://www.mwivet.com).